



## ETHYLENE OXIDE (ETO) GAS CARTRIDGE

We have introduced disposable type 100% ethylene oxide cartridges of various sizes depending on chamber volume. These cartridges are manufactured under stringent conditions. These are available in 40gms, 80gms, 100gms, 170gms, and 300gms capacity of Cartridge.

<b>SR NO.</b>	<b>DESCRIPTION</b>	<b>SIZE</b>	<b>CATRIDGE SIZE</b>
<b>01</b>	<b>EO Aluminum Cartridge (40GMS )</b>	<b>40 gm / 100 Pic</b>	<b>D =35 mm , L =85 mm</b>
<b>02</b>	<b>EO Aluminum Cartridge (60GMS )</b>	<b>60 gm / 100 Pic</b>	<b>D =35 mm , L =85 mm</b>
<b>03</b>	<b>EO Aluminum Cartridge (80 GMS )</b>	<b>80 gm / 100 Pic</b>	<b>D =45 mm , L =150 mm</b>
<b>04</b>	<b>EO Aluminum Cartridge (100 GMS )</b>	<b>100 gm / 100 Pic</b>	<b>D =45 mm , L =150 mm</b>
<b>03.</b>	<b>EO Aluminum Cartridge (170 GMS )</b>	<b>170 gm / 100 Pic</b>	<b>D =53 mm , L =163 mm</b>
<b>04.</b>	<b>EO Aluminum Cartridge (300 GMS )</b>	<b>300 gm / 100 Pic</b>	<b>D =53mm , L =163 mm</b>