

ULTRASONIC CLEANER

ULTRASONIC CLEANER - STF 20

Tank Dimensions : 500 x 300 x 150 mm (LxWxH)

Capacity : 20 Litres (Approx)

MOC : Tank SS304, 16 G and Outer SS

Operating Frequency : 33 KHz \pm 3KHz

Transducer: Imported PZT Sandwich Type Bonded At The Bottom

Ultrasonic Power : 500 Watts

Generator : Digital Microprocessor based Ultrasonic Generator

Digital 0-99 minutes **Timer** with **Temperature display** on LCD

Drain : Provided

Supply : 200-250 V AC , 50 Hz, Single Phase

Accessories : SS basket and lid



ULTRASONIC CLEANER - STF 30

Tank Dimensions : : 500 x 300 x 200 mm (LxWxD)

Capacity : 30 Litres (Approx)

MOC : Tank SS304, 16 G and Outer SS

Operating Frequency : 33 KHz \pm 3KHz

Transducer: Imported PZT Sandwich Type Bonded At The Bottom

Ultrasonic Power : 500 Watts

Generator : Digital Microprocessor based Ultrasonic Generator

Digital 0-99 minutes **Timer** with **Temperature display** on LCD

Drain : Provided

Supply : 200-250 V AC , 50 Hz, Single Phase

Accessories : SS basket and lid

